



Certificate in relation to air transportation of Mobility Aids with Lithium Batteries by passengers with walking disabilities

We hereby certify that our product

e-pilot P15

including the **lithium-ion battery** with article number 1592726 **ID 39533-00** fulfil the requirements of

IATA chapter 2.3.2.4

and therefore the system may be carried as checked baggage without restriction.

Only the battery has to be removed for transportation (protected from unintentionally switching on, damage and from short circuit by other luggage e.g., by placing the battery in a protective pouch) and **carried** in the **passenger cabin**.

The lithium ion batteries are tested according UN Manual of Tests and Criteria" Section 38.3.

The rated energy of the batterie is with 280 Wh **below** the IATA limit of **300 Wh** for batteries which can be removed by the user.

The pilot in command must be informed about the location of the battery.

We strongly recommend contacting the airline sufficiently in advance to request approval for transport.

This certificate does not apply to damaged or defect batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31st December 2024**.

Albstadt, 13th of November 2023

Alber GmbH

A handwritten signature in blue ink, appearing to read "J. Beltman", with a long horizontal stroke extending to the right.

Johannes Beltman
Managing Director